EN

Carefully read the following instructions before use.

Description

All of the BiFlex and BiFlex Evo devices are sterile, for single use and disposable. Please refer to the device labeling for the specifications of each dilator and sheath.

BiFlex and BiFlex Evo are comprised of a flexible, hydrophilic-coated ureteroscopy sheath and a hydrophilic-coated ureteral dilatation catheter with a tapered distal tip and dual working lumen (3 Fr) for guidewires and fluids. This device is available in two models with sheath inner diameters of 10 and 12 Fr.

Indications

The ROCAMED BiFlex and BiFlex Evo are intended to be a conduit for passage of endoscopes and other urological devices for the purpose of performing ureteroscopy procedures. The dual working lumen dilator with luer lock connections allows the user to insert guidewires and fluids.

Contraindications

Presences of a large obstructing distal ureteral stone
Presences of a tight stricture of the ureter

Adverse events

Tissue perforation and tissue lesions.

Precautions

These precautions for use should be fully understood before using the device.

- Do not use the device if damage to the sterile barrier is suspected (i.e. damaged and/or deteriorated packaging).
- Do not use the device if it has reached or passed its use by date.
- Do not use the device if it is damaged upon removal from its packaging.
- Physicians should have thorough understanding of the technical principles, clinical applications and risks associated with ureteral dilatation procedures using sheaths and flexible dilators is required before attempting to use this device.
- The different components of the device and their uses should be properly understood before using the device.
- The sheath and dilator should be handled with the utmost precaution and care should be taken to avoid excessive or unnecessary handling, especially when inserting the assembly and advancing it over the guidewire.
- The device’s hydrophilic coating must be activated before use.
- If resistance is encountered, suspend progression of the device until the cause has been determined.
- If resistance is encountered during withdrawal of the device, suspend withdrawal until cause has been determined.
- In the event of significant stenosis, consider ureteral dilatation prior to insertion of the device.

Directions for use

Remove the BiFlex or BiFlex Evo from its sterile packaging. Prior to insertion of the device, use an endoscope to advance a guidewire, if possible coiled in the renal pelvis. Activate the hydrophilic coating of the device by immersing in sterile physiological saline for approximately thirty seconds. When the guidewire is in place, pass the device over the guidewire. Take care to advance the sheath and dilator jointly; the two devices must be perfectly interconnected (the clip of the dilator must be fitted to the funnel of the sheath). Advance the device slowly so that the dilator can expand the tract as it moves forward. When the device has been positioned and its position confirmed under fluoroscopy, a safety guidewire can be left in place in one of the two working lumens of the dilator. The other lumen is used as needed either for a second guidewire (additional security), to inject contrast medium, or as an access route for a laser fiber or any other endoscopic instrument that can pass through a 3 Fr lumen. To remove the dilator, hold in place the sheath, unclip the dilator from the funnel of the sheath and withdraw the dilator. Once the dilator has been withdrawn, the sheath opening is wide enough to insert a flexible ureteroscope or any other endoscopic devices needed for the procedure. This step can be repeated several times since the sheath thoroughly protects the ureter from damage by the endoscope. Sheath removal is performed by simply pulling gently after having removed any instruments from its working lumen.

Conditions of use and storage

The device should be stored and transported in a normal environment, i.e. away from extreme temperatures and humidity. Do not use if the sterile barrier is damaged and/or deteriorated. These sets are for single use only. Do **not** reuse, reprocess or resterilize. Reuse, reprocessing or sterilization may compromise the integrity of the set and/or lead to failure, which may result in patient injury, illness or death. Also reuse, reprocessing and sterilization may introduce a risk of contamination of the set and/or cause patient infection or cross infection to another patient, which may lead to injury, illness or death. Under any circumstances, the manufacturer will not be held liable in the case of reuse, reprocessing or resterilizing of the device. After use, dispose of set and packaging in accordance with hospital, administrative and/or local government policy.

Warranty and Liability Statement

PROMEPLA and its Affiliates warrant to the first purchaser of this device that reasonable care has been exercised in the design and fabrication of this device. If any suspected damage is found please call the PROMEPLA representative. This limited product warranty is in place of all other warranties, whether express or implied, including but not limited to any implied warranty of merchantability or fitness for a particular purpose. PROMEPLA and its Affiliates under this limited device warranty will be limited to replacement of defective device. **Under no circumstance will ROCAMED Inc., PROMEPLA and its Affiliates be liable for any indirect, incidental or consequential damages resulting from your handling or use of this device.**

		Manufacturer : ROCAMED France – ZI de Signes Allée de Stockholm - 83870 Signes – FRANCE Tel. +33 494 902 100 - Fax. +33 494 986 055
		Distributed by ROCAMED SAM 9 Avenue Albert II, 98000 Monaco Tel: +377 97 98 42 43 Email: info@rocamed.com
		This device is not made with natural rubber latex
Sterilized using ethylene oxide		CE marked on: 20 September 2013

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