EN


Carefully read the following instructions before use.

Description of the device

The Rocamed urethral stent is a sterile single use device comprising two tubular sections in silicone linked by a connecting thread. The stent enables urine drainage outside of the bladder. The two parts of the stent present vertically on each side of the striated sphincter and to avoid migration of the stent. The stent is supplied mounted on its insertion device. The Rocamed urethral stent is for single use (30 days maximum).

EXIME®-M : bladder neck-prostatic apex distance < 5 cm

EXIME®-L : distance apex-coll vesical : 5-7 cm

1. Prostatic part
2. Bulbar part of the stent
3. Connecting thread between the two parts
4. Stabilising wings
5. Safety and removal thread
6. Pusher tube
7. Flattener tube of the wings
8. Flattener tube pink bumper
9. Luer connection between pusher tube and stylet of alignment.

Indications

The Rocamed urethral stent is for drainage of the urine from the bladder via the urethral canal in the adult male in acute or chronic urine retention.

Contraindications

Infection of the genito-urinary system, macroscopic haematuria with clots, sphincter insufficiency. Urinary tract diseases, such as urethral stenosis, urethral false route, bladder stones or other significant diseases that may affect the normal functioning of the device.

Do not use after treatment with a physical agent generating prostatic oedema and possibly macroscopic haematuria (such as, hypo-fractionated radiotherapy, brachytherapy, focused ultrasound, trans-urethral micro-wave thermotherapy...), until the urine becomes clear (48h minimum).

Do not use, when, on the advice of the doctor, such a procedure would be contrary to the greater interest of the patient.

Possible complications

Urinary urgency, weak urinary stream, acute urinary retention, urinary leaks, urethral inflammation, bladder lithiasis, haematuria, acute prostatitis, orchiepididymitis, perineal pain, urinary infection, fever, retrograde ejaculation, and other typical complications linked to the placement of urethral stents.

Adverse events

The adverse events associated with the use of this device include, without being limited to: migration, expulsion, encrustation of the device, breakage of the connection thread between the two tubes of the device, breakage of the removal thread, false urethral route at the time of placement of the device.

Precautions

The following instructions for use must be scrupulously followed in order to ensure correct placement and correct functioning of the device and the safety of the patient:

1. Placement and removal of this device are reserved for practitioners specialised in surgical urology.
2. Do not use the device if the packaging is damaged, compromising the integrity of the latter and/or its sterility.
3. Do not use the device if it has exceeded the expiry date.
4. Check the integrity of the device before use. Do not use in case of suspicion about its integrity.
5. After insertion, regular monitoring of the stent should be put in place by the practitioner to evaluate the performances of the device and check if there are any complications. During this monitoring, any suspicion of abnormality or of incorrect functioning of the device should lead to its removal and implementation of the appropriate treatment.
6. Proceed with caution during insertion and removal of the device. During insertion, stop immediately in the event of excessive resistance to progression of the device in the urethra.
7. Antibiotic prophylaxis is recommended before placement.
8. The use of the device with a magnetic resonance imaging (MRI) device has not been evaluated.
9. Do not use with a patient presenting with allergy or sensitivity to silicone.
10. Do not try to insert a urinary catheter (indwelling Foley catheter) in a patient fitted with this device. In case of acute urine retention, remove the device before insertion of an indwelling Foley catheter.

Instructions for use

The following instructions for use must be scrupulously followed in order to ensure correct functioning of the device and the safety of the patient.

Before use, check the integrity of the device and the correct position of all the elements. Ensure that:

1. The translucent flattener tube covers the wings of the prostatic and bulbar parts of the stent. The pink bumper of the flattener tube must be at the level of the three markers of the pusher tube. The safety and removal thread is correctly blocked in the Luer connection: a slight traction on this thread must show the absence of any mobilisation of the thread.
2. Instill a syringe of lubricating gel and local anaesthetic in the urethral meatus. Stretch the penis vertically and insert the 22 or 24FR bougie carefully so as to dilate the meatus and to check for the absence of urethral stenosis. Insert it until abutment of the tip of the bougie against the posterior wall of the bulbous urethra. Read the insertion on the graduated scale of the pusher tube.
3. After removal of the bougie, keep the penis vertically stretched and insert the device. When the pink bumper of the flattener tube is blocked against the urethral meatus, continue to glide the pusher tube inside the flattener tube while maintaining vertical traction of the penis.
4. Glide the device inside the urethra to the same depth as the one measured previously with the bougie. At this time, the operator feels the abutment of the device against the posterior wall of the bulbous urethra. When the operator feels progression is stopped, the stent is in the correct position. The stylet, the pusher tube and the flattener tube can be removed.
5. Still keeping the penis stretched vertically with the same hand, take hold of the pusher tube in the other hand and withdraw the Luer connector (between thumb and index so as to release the removal thread. Relieve the stylet and the pusher tube vertically.
6. Pull lightly on the removal thread which exceeds the urethral meatus. The perception of an elastic resistance enables to check that the prostatic tube is correctly in place. Its wings abut the striated sphincter of the urethra. If this is not the case, remove the device.
7. Cut the removal thread with a pair of sterile scissors leaving it to come out around 10 mm outside of the meatus.
8. Ask the patient to get up and urinate. The correct position and the correct functioning of the device are confirmed by the absence of leaks in getting down from the examination table or the bed and by the emission of clear urine.

Patient recommendations

It is advisable to explain to the patient using simple and understandable terms the appropriateness of informing his doctor if the patient experiences pain, presents with urinary disorders, fever, irritation of the bladder or any difficulty in urinating.

The practitioner who fits the device should inform the patient of the risks linked to the placement, wearing and removal of the device. He should inform the patient of the procedure to follow in case of anomaly and give the patient the card for users of the EXIME® stent, to keep in his wallet.

Removal:

- It is recommended to calibrate the urethral meatus using a 22FR/24FR bougie before removal to facilitate the passage of the wings of the prostatic tube through the urethral meatus.
- If the removal thread is visible at the meatus: first, take hold of the thread with forceps. Then, instill a syringe of lubricating gel and pull gently on the thread. The perception of resistance is normal at the beginning of removal: it corresponds to the crossing of the striated sphincter by the wings of the prostatic tube.
- If the removal thread is not visible at the urethral meatus: removal must be performed with a fibroscope and forceps for removal of foreign bodies. Progress inside the bulbous urethra and grab either the removal thread or the lower edge of the bulbar tube and remove the stent.
- In case of excessive resistance, stop the removal and determine the cause before continuing.

Conditions of use and storage

Store and transport in a normal environment, i.e., protected from extreme temperatures and humidity. Do not use if the sterility protection is damaged and/or open. Single-use device: the reuse of this device jeopardizes its safety, may generate dysfunction and may also lead to a cross-contamination. In this case, the manufacturer is relieved of any liability. Destroy after use in accordance with the regulation in force on infectious wastes.

Responsibilities and Guarantee

ROCAMED and its affiliates declare to the first-time buyer that reasonable care has been applied in the design and manufacture of this product. In the event of alleged damage, please take contact with ROCAMED and its affiliated companies. This limited warranty acts as any other express or implied warranty, including but not limited to any implied warranty of merchantability or fitness for a particular purpose. The responsibility of ROCAMED and its affiliates under this warranty limited to the product will be limited to the replacement of the product if faulty. On no account may ROCAMED and its affiliates be held responsible for damage, whether direct, indirect and/or incidental resulting from the use and/or the manipulation of this device.

	Manufacturer : ROCAMED France – ZI de Signes Allée de Stockholm – 83870 Signes – FRANCE Tel. +33 494 902 100 Fax. +33 494 986 055	Distributed by ROCAMED S.A.M 9 Avenue Albert II, 98000 Monaco Tel: +377 97 98 42 43 Fax: +377 92 05 61 50 Email: info@rocamed.com
--	---	---

		This device is not made with natural rubber latex	CE marked on: 20 Sept 2013
Sterilized using ethylene oxide		Do not reuse	

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