EN

Description

All ROCAMED RocaFlow Catheters are sterile, for single use and disposable. Please refer to device labeling for the specifications of each model.

Indications

RocaFlow Catheters are indicated for use in the drainage and irrigation of the bladder.

Contraindications

None known

Potential Complications

Irritation of the urethral mucosa, blockage of the catheter due to encrustation, infections. Incorrect positioning of the catheter can cause Urethral damage

Precautions

The following recommendations must be observed strictly in order to achieve maximum device performance and to ensure optimal patient comfort and safety.

- Use only a water-based lubricant
- Inflate the balloon only with Sterile water or Sterile Water with glycerin 9:1 and respect the filling volume mentioned on the balloon. Do not over-inflate.
- Never clamp the Catheter.
- Never use a needle to puncture the catheter shaft for urine sampling.
- Periodic evaluation is advised. Catheter must not remain indwelled more than 2 weeks.
- Once implanted, the physician should monitor the performance of the catheter and check for possible complications. Any suspicion of a defect or incorrect functioning should lead to the removal of the device and onset of appropriate treatment.
- If the balloon loses fluid during implantation, the catheter may move inside the patient. To prevent this risk, the inflation liquid shall be periodically withdrawn and the balloon shall be re-inflated with the stated filling volume.
- Do not use the device if damage to the sterile barrier is suspected (i.e. damaged and/or deteriorated packaging).
- Do not use the device if it has reached or passed its use by date.
- Do not use the device if it is damaged or incomplete either upon removal from its protective sheath or during the procedure.
- Any mishandling or use other than that stated in the instructions for use or contrary to these precautions may result in damage to the device and injuries.
- Physician should have thorough understanding of the technical principles, clinical applications and risks associated with procedures related to this device.
- If resistance is encountered, suspend progression of the device until the cause has been determined.
- If resistance is encountered during withdrawal of the device, suspend withdrawal until cause has been determined.

Instructions for use

Carefully read all instructions prior to using this device. Remove catheter from the peel-open pouch. Lubricate the catheter with a water-based lubricant. Insert the catheter into the urethra. When the desired position is reached, use a Luer syringe to inflate the balloon. Do not over-inflate. Ensure functionality by checking for urine flow from the drainage funnel. Removal: Connect a syringe luer and aspirate gently to remove the inflation fluid from the balloon. Do not apply excessive force on the syringe to avoid a vacuum collapse of the inflation lumen. Deflate the balloon fully before withdrawal.

Conditions of use and storage

The device should be stored and transported in a normal environment, i.e. away from extreme temperatures and humidity. Do not use if the sterile barrier is damaged and/or deteriorated. These devices are for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing of resterilization may compromise the integrity of the device and/or lead to failure, which may result in patient injury, illness or death. Also reuse, reprocessing and resterilization may introduce a risk of contamination of the device and/or cause patient infection or cross infection to another patient, which may lead to injury, illness or death. Under any circumstances, the manufacturer will not be held liable in the case of reuse, reprocessing or resterilizing of the device. After use, dispose of devices and packaging in accordance with hospital, administrative and/or local government policy.

Warranty and Liability Statement

ROCAMED and its Affiliates warrant to the first purchaser of this device that reasonable care has been exercised in the design and fabrication of this device. If any suspected damage is found please call the ROCAMED representative. This limited product warranty is in place of all other warranties, whether express or implies, including but not limited to any implied warranty of merchantability or fitness for a particular purpose. ROCAMED and its Affiliates under this limited device warranty will be limited to replacement of defective device. Under no circumstance will ROCAMED and its Affiliates be liable for any indirect, incidental or consequential damages resulting from your handling or use of this device

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