Carefully read the following instructions before use.

Description of the JFil and MiniJFil stent sets

The Rocamed JFil and MiniJFil ureteral stents are sterile and single-use. Refer to the packaging for the specifications of each product.

JFil & MiniJFil: Ureteral stent with a single renal loop and a 0.3 Fr surgical thread for maximum 365 days of continuous use.

Indications

The Rocamed JFil and MiniJFil ureteral stents are used for temporary internal drainage of the ureteropelvic junction or for the preparation of the ureter before a medical intervention. The device has to be implanted at least 15 days before intervention. The length of the product is appropriate if the stent's tail extends at least 3 cm below the obstacle. Use the MiniJFil only in a fully healthy ureter.

The stents may be placed using endoscopic techniques.

The stents are recommended as permanent indwelling devices, it is recommended that the indwelling time does not exceed 365 days.

Target population: adults only (at least 22 years old)

Contra-indications

Untreated active infection of the higher tract. Ureteral perforation. Significant bleeding of the urinary tract.

Immune chronic disease. Impossibility of using NSAIDs or corticoids.

Do not attempt stent placement in a patient with suspect ureteral avulsion.

Do not use, when, in the judgement of the physician, such a procedure would be contrary to the best interest of the patient.

Possible complications

Conventional complications related to the placement of ureteral stents and to urologic procedures.

Adverse events

Complications associated with the use of ureteral stent should include but are not limited to: migration or dislodgement, formation of knots, encrustation, infection, fragmentation, flank pain, mucosal irritation or inflammation, frequency, urgency, dysuria, stone formation, obstruction, erosion.

Precautions for use

The following precautions must be taken scrupulously in order to ensure the proper functioning of the device and the safety of the patient:

1. The use of this device is limited to practitioners trained in urologic surgery procedures and in the placement of ureteral stents.
2. Do not use the device if the packaging is damaged, thus compromising the integrity of the latter and/or its sterility.
3. Do not use the device if it has exceeded the limit date of use.
4. Before implantation, make sure that the dimensions and the type of the selected stent are adapted to the patient (see indications).
5. Control the condition of the device before any use. Do not use in case of doubt about its integrity.
6. Control during implantation time: after implantation, a regular monitoring of the stent must be carried out by the practitioner in order to assess its performance and check any potential complications. During the monitoring, any suspicion of a fault or malfunction of the device must lead to its removal and the establishment of a new treatment.
7. The user must make sure that the dimension of the stent is adapted to the patient. In particular, the device is not intended for pediatric use.
8. Insert the device with caution. Stop immediately in case of anomaly or excessive resistance. Determine and eliminate the causes of the anomaly before continuing.
9. Use only with the provided pusher.

NOTA BENE: The guidelines provided in the sets are those recommended by the manufacturer for the placement of the devices. The use of other guidelines could compromise the procedure. For further information, refer to the Rocowire guidelines for use.

For the models delivered without a guidewire, the manufacturer recommends the use of guidewires that are equivalent to the Rocamed models (hydrophilic wires, stent, with a maximum diameter of 0.035" and a length of 1.5m).

The user must check the compatibility of the guidewire and the stent before use.

Instructions for use

The following precautions must be taken scrupulously in order to ensure the proper functioning of the device and the safety of the patient:

For models with hydrophilic treatment: activate the treatment by immersing the stent in sterile water or physiological saline solution for approximately 30 seconds.

Models with connectable pusher:

1. Insert the ureteral stent over the guidewire, with the tapered tip first.
2. Put the pusher supplied with the stent over the guidewire and insert its end into the beveled tail of the stent. Do not exert excessive force on the stent or pusher. The stent should approximately be inserted up to the black mark printed on the pusher.
3. Move the pusher over the guidewire in order to push the stent up to the desired place. Check the correct position of the two ends by X-ray or endoscopic control.
4. Slowly remove the guidewire while holding the pusher in position and check that the distal loop closes in the desired place before proceeding with the complete withdrawal of the guidewire.
5. Unscrew the two parts of the pusher and separate them in order to release the stent. In the event of resistance, do not force.
6. Remove the pusher from the patient.
7. Perform a final control to ensure that the stent is placed correctly, and the loop is closed properly. Check that the surgical thread emerges through the ureteral meatus. For women, slightly pull on the thread to make sure it is tight, then cut it flush with the urethral meatus.

Models with simple pusher:

1. Insert the stent over the guidewire, orienting the side with the metal ring toward the other side in contact with it.
2. Move the pusher over the guidewire in order to push the stent up to the desired place. Check the surgical thread remains unfolded and knot-free when pushing up the stent. Confirm correct position of the two ends by X-ray or endoscopic control.
3. Confirm correct position of the two ends by X-ray or endoscopic control. Slowly remove the guidewire while holding the pusher in position and check that the distal loop closes in the desired place.
4. Remove the pusher from the patient.
5. Perform a final control to ensure that the stent is placed correctly and the loop is closed properly. Check that the surgical thread emerges through the ureteral meatus. For women, slightly pull on the thread to make sure it is tight, then cut it flush with the urethral meatus.

Withdrawal

Retire cistoscopically by gently pulling on the surgical thread, and if resistance is encountered during removal, fluoroscopically determine stent position and cause of resistance.

Remove the stent with an appropriate grasper. Do not use an instrument that could damage the surgical thread or the implant. If any resistance appears, halt the withdrawal of the device and assess the causes before continuing. If necessary, continue the withdrawal with another method.

Recommendations for patients

The practitioner must inform the patient about the risks associated with the installation of a JFil or MiniJFil ureteral stent, and that the practice of violent sports is contraindicated during the duration of the implantation. He must indicate to the patient the procedure to follow in case of anomaly.

The withdrawal of the stent must be performed according to the instructions of a practitioner and under his supervision.

Conditions of use and storage

Store and transport in a normal environment, i.e. protected from extreme temperatures and humidity. Do not use if the sterility protection is damaged and/or open. Single-use device: the reuse of this device jeopardizes its sterility, may generate dysfunctions and may also lead to a cross-contamination. In this case, the manufacturer is relieved of any liability. Destroy after use in accordance with the regulation in force on infectious wastes.

Responsibilities and Guarantees

ROCAMED and its affiliates declare to the first-time buyer that reasonable care has been applied in the design and manufacture of this product. In the event of alleged damage, please take contact with ROCAMED and its affiliated companies. This limited warranty acts as any other express or implied warranty, including but not limited to any implied warranty of merchantability or fitness for a particular purpose. The responsibility of ROCAMED and its affiliates under this warranty limited to the product will be limited to the replacement of the product if faulty. On no account may ROCAMED and its affiliates be held responsible for damage, whether direct, indirect and/or incidental resulting from the use and/or the manipulation of this device.

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