EN

Carefully read the following instructions before using the RocaNepbro Needles.

Description

All ROCAMED RocaNepbro Needles are sterile, for single use and disposable. Please refer to device labeling for the specifications of each model.

Indications

RocaNepbro Needles are intended for initial puncture during percutaneous nephrostomy procedures.

Contraindications

Major homeostasis issues

Potential Complications

Infection, internal bleeding and / bleeding at the puncture site, puncture or injury to target organ or to a nearby organ, peritonitis, damage of tissue adjacent to the needle.

Precautions

The following recommendations must be observed strictly in order to achieve maximum device performance and to ensure optimal patient comfort and safety.

- Do not use the device if damage to the sterile barrier is suspected (i.e. damaged and/or deteriorated packaging).
- Do not use the device if it has reached or passed its use by date.
- Do not use the device if it is damaged or incomplete either upon removal from its protective sheath or during the procedure.
- Handle this device carefully to avoid damages to any by-stander.
- Any mishandling or use other than that stated in the instructions for use or contrary to these precautions may result in damage to the device and injuries.
- Precautions must be observed to prevent damages to guidewires if advanced or withdrawn through this device. This device is not adapted to guidewires bigger than 0,035".
- Needles should be handled and used with care, under ultrasonic or radioscopic monitoring to prevent damage to the needle and to ensure perfect atraumatic progression towards the target area. If there is the slightest doubt, the physician should withdraw the device under ultrasonic or radioscopic guidance with as much care as possible.
- This needle is not suitable for use under MRI (Magnetic Resonance Imaging)
- Needles should be handled in aseptic conditions, in an operating room, according to the regulations in force in each country.

Instructions for use

Carefully read all instructions prior to using these needles.

Remove the RocaNepbro needles from the peel-open pouch. Remove the needle protective cap and set the sliding stopper to the desired penetration distance. Introduce the needle to the desired distance. During introduction, hold firmly the handle of the cannula and the hub of the stylet to prevent disconnection. Remove the stylet to ensure that desired area has been reached. Remove the inner cannula to allow introduction of a guidewire. Withdraw the needle and dispose in a container in accordance with local regulations.

Conditions of use and storage

The device should be stored and transported in a normal environment, i.e. away from extreme temperatures and humidity. Do not use if the sterile barrier is damaged and/or deteriorated. These devices are for single use only. Do **not** reuse, reprocess or resterilize. Reuse, reprocessing of resterilization may compromise the integrity of the device and/or lead to failure, which may result in patient injury, illness or death. Also reuse, reprocessing and resterilization may introduce a risk of contamination of the device and/or cause patient infection or cross infection to another patient, which may lead to injury, illness or death. Under any circumstances, the manufacturer will not be held liable in the case of reuse, reprocessing or resterilizing of the device. After use, dispose of devices and packaging in accordance with hospital, administrative and/or local government policy.

Warranty and Liability Statement

ROCAMED and its affiliates declare to the first purchaser of this device that reasonable care has been taken in the design and manufacture thereof. If damage is suspected, please contact ROCAMED and its representatives. The present limited guarantee dispenses all other guarantees, express or implied, including but not limited to any implied guarantees concerning merchantability or suitability for specific purposes. The liability of ROCAMED and its companies affiliated by virtue of this limited product guarantee shall be limited to the replacement of defective products. Under no circumstances shall **ROCAMED and its affiliates be held liable for direct, indirect, incidental and/or consequential damages of any kind arising from the use and/or handling of this device.**

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