**Description**

The irrigation/aspiration lines are sterile and for single use only. Please refer to packaging for the specifications of each product.

**Indications**

Single-use sterile system intended for urology, laparoscopy, hysteroscopy, or arthroscopy procedures for insufflation (creation of pneumoperitoneum in laparoscopy), irrigation/washing and aspiration of fluids in anatomical cavities. Refer to the product packaging and/or technical data sheet.

**Contraindications**

These lines should be used after an assessment by the surgeon and under his or her own responsibility. These sets are not suitable for the transfer of blood or blood derivatives. These sets must be used only with fluids or medications approved for the surgical procedure.

**Adverse effects**

Risk of infections. Risk of air embolism following the CO<sub>2</sub> injection during laparoscopic surgery. Loss of visibility of anatomical structures due to excessive use of fluids/medications.

**Precautions for use**

The operator should observe the following precautions conscientiously before using the set:

- 1. Do not use the set if the packaging is damaged because this could compromise device integrity and/or sterility.
- 2. Do not use the set if it has reached or passed its expiry date.
- 3. Do not use the set if damaged when removed from its packaging.
- 4. A thorough understanding of the necessary principles and techniques relating to the lines' use as well as the clinical applications and associated risks is required before use.
- 5. A thorough understanding of the various components of the lines and their use is required before using this set.
- 6. Inspection of the CO<sub>2</sub> injection during laparoscopic surgery.
- 7. For any use with Endoflow® machines, refer to their IFU.

**Instructions for use**

Please open the packaging following aseptic procedures

A. For lines used with the Endoflow® or Endoflow® II pump:
 

- Place the bag(s) in the Endoflow® or Endoflow® II pump;
- Spike the bag(s) port
- Ensure there are no kinks in the tubing
- Close the device doors
- Select the appropriate pressure and/or program if necessary
- Prime the line and check correct operation of the system before use on the patient and outside the patient

B. For lines used under gravity:
 

- Place the bag(s) on the drip stand at the required height
- Spike the bag(s) port
- Ensure there are no kinks in the tubing
- Prime the line and check correct operation of the system before use on the patient and outside the patient

If the tubing is coloured, the blue line is always for irrigation, and the red line is always for aspiration. If a handle is present, particularly for laparoscopy procedures, the blue button controls the irrigation, and the red button controls the aspiration.

**C. Recommendations for use of the aspiration function:**

- To prevent the backflow of liquid, it is recommended to fully drain fluids from the cannula before releasing the aspiration button (about 2 or 3 seconds)
- In order to prevent the aspiration from getting clogged and/or the aspiration piston from getting blocked, it is recommended not to aspirate objects which diameter is greater than that of the cannula.
- During a long operation or when the aspirated fluids are particularly polluted, it is recommended to perform a flush by pressing both pistons simultaneously for several seconds.
- To ensure satisfactory aspiration, it is necessary to verify that the operating room's aspiration system meets accepted standards.

**Conditions of use and storage:**  
The set should be stored and transported in a normal environment, i.e. away from extreme temperatures and humidity. Do not use if the sterile packaging is damaged and/or deteriorated. Single-use system: do not reuse, reprocess or resterilize. This could compromise integrity of the system and/or cause malfunctions potentially leading to injury, illness or death of the patient. In addition, there is a risk of contaminating the system, which could cause an infection in the patient or cross-infection with another patient which could result in injury, illness or death. In this case, the manufacturer will not be held liable. After use, discard the device and packaging in compliance with hospital, government and/or local authority policies.

**Liability and warranty:**  
PROMEPLA and affiliates declare to the first purchaser of the device that reasonable care was taken in the design and manufacture of this device. If damage is suspected, please contact PROMEPLA and its representatives. The present limited warranty dispenses all other warranties, express or implied, including but not limited to any implied warranties concerning merchantability or suitability for specific purposes. This limited guarantee ensures that the responsibility of PROMEPLA and affiliated companies will be limited to the replacement of defective products. **Under no circumstances shall PROMEPLA and affiliates be held liable for direct, indirect, incidental and/or consequential damages of any kind arising from the use and/or handling of this device.**

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