## ROCAMED =

# **KOBOT**Nitinol Stone Retrieval Tipless Basket

- A new range of Stone Retrieval Baskets in line with the concept Of Cobotic. These baskets are diversified, effective and integral to the working of the surgeon and his assistants and allows for maximize performance.
- An innovative and ergonomic Criterium Handle with assisted closing system, helps in catching, retaining and removing of stone fragments.



\*> \*> \*> ENDO-UROLOGY



### **KOBOT**™

## Nitinol Stone Retrieval Tipless Basket





#### **KOBOT TIPLESS**

4 wire basket
 Assisted Closing



#### **KOBOT SIX**

• 6 wire basket



#### **KOBOT FILTER**

• 4/16 wire filter



#### **KOBOT MITT**

• 3 wire



Available with assisted closing System





#### **→** FLEXI TIP IMPROVEMENT

For a better scope deflection



#### **Ordering Informations**

Product code	Size Fr	Wires	Working Lenght (cm)	Basket Ø (mm)	Basket Lenght (mm)	Description
ROSZ1930ST	1.9	4	120	12	17.5	Kobot Tipless with criterium handle
ROSZ1940ST	1.9	4	120	12	17.5	Kobot Tipless - Criterium handle with assisted closing system
ROKT1800ST	1.8	4/16	120	14	18	Kobot <b>Filter</b> with criterium handle
ROKS1700ST	1.7	6	120	10	15	Kobot <b>Six</b> with criterium handle
ROKM1900ST	1.9	3	120	10	10	Kobot MITT with criterium handle

MKG-FY-KOBOT-EN-03
Designed by FOCALYS. Illustrations in this document are non-contractual.
ROCAMED reserves the right to change product specifications without notice.

For more information visit www.rocamed.com







