



ROCAMED 





RIO TRACER™



» » » ENDO-UROLOGY

www.rocamed.com





➔ Hydrophilic Floppy Tip

Facilitates passage beyond obstructions and negotiates tortuous anatomy

➔ STIFF Nitinol Core

Offers a kink resistant shaft for enhanced control

➔ PTFE Coating


Offers a wire surface to facilitate advancement

➔ Blue and Yellow striped pattern

Designed to provide clear endoscopic visualization of wire movement

➔ Barium Filled Radiopaque Tip

Enhances Fluoroscopic Visualization



Ordering Informations

Product code	Ø (inches)	Length (cm)	Distal Tip (3cm)
ROWR3502ST	Ø 0,035	150	Hydrophilic and Floppy

MKG-FY-MEN-RIOTRACER-EN-02
 Designed by FOCALYS. Illustrations in this document are non-contractual.
 ROCAMED reserves the right to change product specifications without notice.

For more information visit www.rocamed.com