

# SIRIUS

60W THULIUM FIBER LASER




## STONE & BPH

- EXTREME FREQUENCY EMISSION
- LOW RETROPULSION
- EFFECTIVE LITHOTRIPSY
- GREAT CUTTING & ABLATION CAPABILITY






# SIRIUS



## STONE SOLUTION

## BPH SOLUTION



Reduced Size of Dust



Low Retropulsion



Extreme Frequency Emission



Compact Design



Reduced noise levels



Standart Electrical Installation



BPH



Hemostasis




Soft Tissues




## FIBERS

Single Use & Reusable



## FOOT PEDALS



## INTUITIVE SCREEN





# SIRIUS

## SOLUTION




### Technical Informations


Laser type	Thulium Fiber Laser
Wavelength	1940 nm $\pm$ 20nm
Power of Target	60 Watts
Energy/Pulse	Up to 6J
Application mode	Pulsed and CW
Time function	Pulsed Mode and CW
Pulse Duration	Up to 15ms
Repetition rate	Up to 2500 Hz
Peak power	Up to 500W

### Ordering Informations

Reference	Description
MS01	SIRIUS 60w Thulium Medical Laser
FIBERS	Single Use & Reusable

For more information visit [www.rocamed.com](http://www.rocamed.com)

 Rocamed SAM  
9 avenue Albert II  
98000 Monaco

 (+377) 97 98 42 43

 [info@rocamed.com](mailto:info@rocamed.com)

 [www.rocamed.com](http://www.rocamed.com)