

TraxerFlow[™] II Irrigation Line



→> →> UROLOGY

www.rocamed.com

TraxerFlow[™] Irrigation Line

ROCAMED

TraxerFlow[™] Dual Port

→> PRESSURE CHAMBER

• Anti-reflux valve prevent back flow of fluid.

MANUAL PUMP

• Provide extra flow for clearing the surgical field.

➔ LUER LOCK CONNECTION

• Grips the devices and prevent leakage.

→> REMOVABLE TUBING EXTENSION

• Allow surgeon to use the pump itself or with the assistant.

→> TRANSPARENT BLUE TUBE

• Allow to see the saline water.

Product description

GRAVITY LINE

Référence

ONGT2000ST12 TraxerFlow DUAL PORT Irrigation Line for Ureteroscopy Laser







Compatible with all principal brands of endoscopes

WITH ENDOFLOW II

Référence	Product description	
SOUR9107ST12	TraxerFlow Irrigation Line for Ureteroscopy Laser	1 spike
SOUR9207ST12	TraxerFlow Irrigation Line for Ureteroscopy Laser	2 spikes















MKG-FY-TFLOW-EN-00 Designed by FOCALYS. Illustrations in this document are non-contractual. ROCAMED reserves the right to change product specifications without notice.

For more information visit www.rocamed.com



Rocamed SAM 9 avenue Albert II 98000 Monaco







