


ROCAMED 

EndoFlow™ II in Gynecology

Hysteroscopy and Laparoscopy



⇒ User friendly manipulation

⇒ Self-regulation and accurate pressure management



⇒ ⇒ ⇒ GYNECOLOGY

www.rocamed.com



EndoFlow™ II

NEW GENERATION,
A UNIQUE SYSTEM OF FLUID MANAGEMENT


- » ACCURATE MANAGEMENT OF THE PRESSURE
 - Self-regulation
 - No pulsation effect
- » SIMPLICITY AND USER FRIENDLY MANIPULATION
- » PRESSURE UNDER CONTROL
 - Tubings RFID managed
- » BAG HEATED AND MAINTAINED AT 38° ±2



HIGH PERFORMANCE IRRIGATION




HANDPIECE:

- Red piston for suction control
- Free pressure regulation
- Left and right hand use



PROBES:

- Straight liquid jet
- Insulating treatment
- DEHP free / Phtalate free



LAPAROSCOPY TUBING SET (IRRIGATION/ASPIRATION)

References	Description	N. Spike
SOLA0533ST12	Handpiece + Probe (for Ø 5mm probe tip; L:330mm)	1 spike
SOLA1033ST12	Handpiece + Probe (for Ø 10mm probe tip, L:330mm)	1 spike

HYSTEROSCOPY TUBING SET

References	Description	N. Spike
SOHY9102ST12	Aspiration tubing (push fitting)	1 spike

MKG-FY-MEN-GYNECO-EN-02
Designed by FOCAYS. Illustrations in this document are non-contractual.
ROCAMED reserves the right to change product specifications without notice.

For more information visit www.rocamed.com